3.2.S.2.4. HOLD TIMES [BNT MAINZ AND RENTSCHLER]

The validated hold times for manufacturing are described in Table 3.2.S.2.4-1. Pool hold times are not required for routine processing, but strategic holds in the process (≥24 h) to aid in manufacturing scheduling were validated as described in Section 3.2.S.2.5 Hold Times [BNT Mainz and Rentschler].

Table 3.2.S.2.4-1. Hold Times for mRNA Drug Substance (Routine Manufacturing)

4.2 1st ind

During the Proteinase K Pool hold time the transport at 4.2 1st ind to Rentschler, Laupheim for further processing is included (Section 3.2.S.2.5. Shipping Performance Qualification [BNT Mainz and Rentschler]).